Allium Bulbar Urethral Stent (BUS) System

- **Bulbärer Harnröhren-Stent (BUS) System von Allium- DE**
- **Allium système de Stent Urétral Bulbaire BUS- FR**
- **Sistema Stent da Uretra Bulbare (BUS) Allium- IT**
- **El Sistema de Stent Uretral Bulbar (BUS) de Allium- SP**

- Instructions for Use -

- **Bedienungsanweisungen –DE**
- **Mode d'emploi-FR**
- **Istruzioni per l'Uso-IT**
- **Instrucciones de Uso-SP**

Manufactured by

Allium Ltd.

Allium Bulbar Urethral Stent (BUS) System

Device description:

The Allium **Bulbar Urethral Stent (BUS)** is a temporary implant intended for trans-urethral insertion into the male bulbar urethra diagnosed with a stricture in order to open the occluded passage and allow spontaneous urination. The BUS is made of a super-elastic coil covered with a co-polymer. Once inserted into the bulbar urethra by a delivery system, the stent is released to allow its self-expansion. The stent is intended to remain in place up to 3 years.

Table 1 – BUS Dimensions

Nominal BUS stent diameter (Fr)	Nominal BUS stent lengths (mm)	Profile of delivery system (Fr)	Working Length of delivery system (cm)
(45)	50/60/70/80	22	20

The BUS set consists of:

1. Bulbar Urethra Stent (Figure 1)
2. An Endoscopic Delivery tool (Figure 2).
3. Meatal Shield (Figure 3).

Figure 1 – Bulbar Urethra Stent (BUS)

Figure 2 - Endoscopic Delivery Tool



Figure 3 - Components of the Meatal shield

- a - Assembled Meatal Shield
b - Dilator Sheath
c - Dilator shaft \ mandrel

The Meatal Shield minimizes trauma to the distal urethra and meatus during insertion of the delivery system. The Meatal shield (a) is composed of a Dilator Sheath (b) and a Dilator shaft (c).

**Indications for Use:**

Management of Bladder Outlet Obstructions (BOO) caused by bulbar urethral strictures.

Contraindications:

The BUS is contra-indicated in the following cases:

1. Are younger than 18.
2. Have an acute infection of the prostate or urethra.
3. Febrile patients
4. Active urinary tract infection
5. Immuno-compromised patients, who have a prosthetic heart valve or other implanted device
6. Patients who cannot tolerate any form of antibiotic treatment
7. Irreversible atonic/acontractile bladder
8. Anticoagulation therapy and/or bleeding disorders
9. Urethro-cutaneous/urethro rectal fistula
10. Present urinary calculi or Cystolithiasis
11. "Constriction rings" and/or vacuum erection devices
12. Use injectable medications to obtain erection
13. Patients with penile implants and/or implanted artificial urinary sphincters
14. A history of any illness, medications or surgery that may affect the efficacy of the stent

Informing Patients, Potential Adverse Events and Complications:

The following adverse events might occur due to the implanted stent or the implantation and extraction procedures:

1. Following the procedure some patient may endure post-voiding urinary dribbling and post-ejaculatory seminal dribbling.
2. Some patients can develop reactive urinary frequency or urgency
3. Reactive tissue proliferation at one or both ends of the stent. In most cases, this tissue disappears spontaneously after stent removal.
4. Hematuria
5. Pain/discomfort
6. Bleeding
7. Urinary tract infection
8. Stent migration. If this occurs, the stent should be removed and replaced with a new one
9. Stent encrustation. If this occurs, a lithotripsy should be performed or the stent should be removed
10. Stent misplacement
11. Mechanical damage to the urinary tract

Precautions:


1. Do not use after the “use by” date.
2. The packaging should be inspected for damage.
3. The delivery system should be inspected for damage prior to use.
4. Prior to delivery, the physician should ensure that the stent is covered by the over-tube.
5. The BUS system should be used only by physicians trained in interventional or endoscopic techniques.
6. Prior to its use, the technical information supplied with the device should be carefully reviewed.

Warnings:

1.  The BUS is intended for **Single Use Only-DO NOT RE-STERILIZE**. Reuse, reprocessing, re-sterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics and may lead to device failure, which may result in injury to the patient. It may also create a risk of contamination of the device which may lead to injury, illness or death of the patient or end user.
2. The product should be carefully inspected prior to use to verify that the sterile package has not been damaged. The device should not be used if the package is open or damaged, or if the device has been contaminated prior to insertion.
3. Remounting the stent onto the delivery system should not be attempted.
4. Positioning of the main body of the stent across the external sphincter may render the patient incontinent.
5. The stent is not intended as a permanent option to treat urethral obstructions.
6. Bladder catheterization with an implanted stent is not recommended. Introduction and passage of a catheter through the stented urethra into the bladder may dislodge the stent and/or damage the polymeric cover.
7. The use of transurethral instrumentation while the stent is in place is not recommended because longitudinal compression of the stent by instrumentation could dislodge the stent.
8. The stent should be inserted and positioned only according to the Instructions for Use (IFU) that are detailed at this document.
9. Under-vision insertion should not be attempted in patients in whom bleeding may impede the visualization process.
10. After use, the BUS may be a potential biohazard. Handle and dispose of it, in accordance with any local and federal laws regarding medical waste.

PROCEDURE

Figure 4 – Allium BUS segments and correlation with the bulbar urethra anatomy

- * The "body" of the stent is composed of the high radial force "Main Body" and the "Sphincteric Segment" having a gradual decreasing radial force. The main body diameter is 45 Fr, the "Sphincteric segment" has decreasing diameter down to 30 Fr.
- ** The BUS lengths detailed above comprise of the "Main Body" and the Sphincteric Segment combined.

Pre Treatment Evaluation of Bladder Outlet Obstruction (BOO):

It is recommended to perform the following basic evaluation prior to initiation of treatment for BOO with a stent: history, physical examination, PVR (post-void residual urine), DRE (digital rectal examination), prostate sonography (abdominal and/or trans-rectal), uroflowmetry, PSA (prostate specific antigen), urethrography and urinalysis.

If indicated by the results of the basic evaluation, a more extensive workup may be required.

Patient preparation

Prescribe a broad-spectrum oral antibiotic treatment before the insertion and continue this treatment in combination with a non-steroid anti-inflammatory drug according to hospitals endoscopic transurethral procedures protocol.

Equipment required in the operating room:

1. 4 mm (12 Fr) rigid Optic Element (0°, 5° or 30°) with working length of 280 to 305 mm.
2. X-ray equipment
3. Dilatation kit or equipment for urethrotomy
4. Lubricating gel
5. Local anesthetic gel
6. Measuring devices

Measurement of the stenosis

1. Measure the stricture length from its distal end to the sphincter and add 1 cm, the Sphincteric Segment of the Stent will be situated just below the external sphincter and the "Body" segment of the stent should cross the proximal end.
2. Choose the appropriate stent length (50/60/70/80mm).

Dilation of the stenosis

The stenotic area must be dilated to **28 Fr** or incised (Y-shape incision).

Preparation of the BUS delivery system


1. Remove the delivery system from its tray.
2. Insert the Cystoscope optical element into the delivery system until the lens reaches the edge of the delivery system. Optimal optical positioning is achieved when a crescent shape is visible at the bottom of the optic view (see Figure 5).

Figure 5 – Optimal optical element positioning


3. The optical element is held in place. In case that there is a need, repositioning of the optical element can be obtained simply by moving it along the optical element channel until the crescent shape is visible again.
4. Connect the irrigation fluid set to the inlet port (see A in Figure 6).
5. Connect the light cable to the light cable port of the Cystoscope (see B Figure 6).

Figure 6 – Cystoscope optical element insertion

280 mm Cystoscope

300 mm Cystoscope

Inserting the Meatal Shield

1. Check that the mandrel of the Meatal Shield is fully inserted into the Sheath.
2. Lubricate the Meatal Shield and gently insert it into the urethral meatus in a rotating motion.
3. Retrieve the mandrel leaving the Dilator Sheath in position (see Figure 3 (b) and Figure 3 (c)).

Stent insertion & positioning

1. Insert the delivery system with the optical element through the Meatal Shield tube and advance it under vision until you see the distal end of the external sphincter.
2. Assure that the delivery system does not penetrate the sphincter and that the stent covers the entire stricture


Stent deployment

1. Unlock the trigger by pressing the safety lock button (see Figure 2).
2. While keeping the delivery system in the desired location, start squeezing the trigger gently for stent deployment. Pull the trigger all the way repeatedly at least 10 times depending on stent length. The Teflon sheath should gradually come out from the posterior end of the delivery tool.
3. In case that the Teflon sheath does not come out gradually from the posterior end and/or pulling the trigger becomes difficult after a few trigger pulls, please follow the instructions below:
 - Maintain a fixed position at all times
 - Grip tightly the Teflon Sheath button located at the posterior side of the handle (see Figure 7).
 - Pull the button backwards.
4. Wait for 2 minutes for the stent to expand.
5. Acquire a fluoroscopy/X-ray image and make sure the stent is expended, in place and is fully released from the delivery system.
6. Under vision, slowly pull the delivery system outward with semi-circular movements while verifying that the stent maintains its position.
7. Verify that the stent is not penetrating the sphincter (see figure 8).

Figure 7

Figure 8 - BUS accurate position

Removing the Stent

1. The stent can be removed under vision using a rigid endoscopic or a strong flexible endoscopic foreign body forceps.
2. Under vision, insert a 19-21 Fr Cystoscope until the down-stream end of the stent is seen in the bulbar urethra.
3. Insert the foreign body forceps through the Cystoscope.
4. Engage one of the stent distal metal edges (see Figure 9) using the forceps.
5. Start pulling the metal edge outwards. This may initiate tearing of the polymeric cover. Continue to pull outward the metal wire together with the endoscope sheath.
6. Verify that the entire stent came out by checking that the loops at both stent edges were retrieved (see Figure 9).

Figure 9 – BUS stent**MRI Information**

MR Conditional

The Bulbar Urethral Stent was determined to be MR-conditional.

Non-clinical testing demonstrated that the Bulbar Urethral Stent is MR Conditional. A patient with this implant can be scanned safely, immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-Related Heating

In non-clinical testing, the Bulbar Urethral Stent produced the following temperature rises during MRI performed for 15-min of scanning (i.e., per pulse sequence) in 1.5-Tesla/64-MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems:

	1.5-Tesla	3-Tesla
MR system reported, whole body averaged SAR	2.9-W/kg	2.9-W/kg
Calorimetric measured values, whole body averaged SAR	2.1-W/kg	2.7-W/kg
Highest temperature change	+2.0°C	+2.3°C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Bulbar Urethral Stent. Therefore, optimization of MR imaging parameters to compensate for the presence of BUS stent may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 5-mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	1,212-mm ²	231-mm ²	1442-mm ²	317-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

DISCLAIMER OF WARRANTIES

Allium, Ltd. warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including but not limited to any warranties of merchantability or fitness for a particular purpose. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact, and since Allium, Ltd. has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after the device leaves our possession, Allium, Ltd. does not warrant either a good effect or against any ill effect following its use. Allium, Ltd. shall not be liable for any incidental or consequential loss, damage, or expense arising directly or indirectly from the use of this device. Allium, Ltd. will replace any device that we feel was defective at the time of shipment. No representative of Allium, Ltd. may change any of the foregoing or assume any additional liability or responsibility with this device.

Labeling Information

Symbol	This Symbol Means
	Single Use
	Use By
	Batch Code
	Sterilization Using Ethylene Oxide
	Catalog Number
	Caution, Consult Accompanying Documents
	Manufacturer
	Authorized Representative in the European Committee
	Consult Instructions for Use
	Store in Dry Place at Room Temperature
	Do Not Use if Package is Damaged

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MEDNET GmbH