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DESCRIPTIONS

FIAB vessel dilators are intended for single use and are supplied sterile and non-pyrogenic.

INDICATIONS

The FIAB vessel dilators are designed for percutaneous entry into a vessel in order to enlarge the opening of the vessel for the placement of a catheter in a vein or artery.

CONTRAINDICATIONS

The product may be contraindicated in case of known or suspected vascular obstruction. Relative contraindications may include patients with prior irradiation of the subclavian vein, severe scoliosis and chest traumas.

WARNINGS

- FIAB dilators are designed for one time use only. **DO NOT RE-USE**. Reusing the product may lead to: possible presence of biological residues that may cause cross-infections; alteration of materials; loss of initial functional features of product.
- Do not use if package is opened or damaged.
- The product is sterile and non-pyrogenic in unopened and undamaged package. **STERILIZED BY ETHYLENE OXIDE. DO NOT RESTERILIZE**. Resterilization could change the physical characteristics of the material and should not be attempted. Do not sterilize in an autoclave.
- Do not use the device after the expiration date written on the package.
- Caution when using this device. Be aware of needles and sharp objects used during the procedure.
- Do not use the device if it is damaged or deformed.
- Do not overtighten. Do not proceed if resistance is felt or if the components happen to disconnect.

FIAB vessel dilators are designed for use with guidewires. The maximum diameter of guidewire to be used is specified on individual vessel dilator labels. The use of guidewires has been associated with greater incidence of thrombus. Optimal guidewire size and judicious use are recommended.

Variations in individual patient anatomy may preclude the utilization of the percutaneous technique.

PRECAUTIONS

- The device must be used only by or under the supervision of a physician trained in percutaneous vascular access procedures.
- Employ aseptic technique during removal from the package and during use.
- To avoid damage or crimping of the dilators during removal from the package, grasp the fitting element and slowly withdraw from the package.
- Care should be taken during placement and withdrawal of dilator in order to prevent possible tissue damage.
- Pay attention to the guidewire being used during the procedure: for the 4.0FR size FIAB dilator the correct guidewire size to be employed is the 0.025" one; for dilators of different size, the 0.038" is the appropriate choice. FIAB takes no responsibility for any misuse of the dilator.

DIRECTIONS FOR USE

For the procedure to use please refer to the Seldinger technique. The following instructions are provided as a general guide intended for informational purposes only; physicians should add to or alter procedural details by considering and following their own clinical experience and expertise.

- 1) Identify the correct blood vessel
- 2) Insert needle into the chosen vessel according to the standard percutaneous technique.
- 3) While firmly holding needle in place, insert guidewire through the needle cannula until desired spot is reached.
- 4) Once the position of guidewire is verified by fluoroscopy or appropriate technique, pull the needle out while holding the wire in place.
- 5) Place dilator over guidewire to dilate the blood vessel tissue.
- 6) To insert the desired device or catheter, follow manufacturer's instructions.

STORAGE

The product must be stored in its original package at the environmental conditions (temperature and relative humidity) specified on the pouch's label. Putting external heavy weights on the package could damage the product.

WASTE DISPOSAL

Waste coming from health facilities must be disposed of according to the regulation in force.

GENERAL NOTES

If a serious incident occurs while using this device or as a result of its use, please report it to the manufacturer and to your national authority. For any malfunction or defect of the device, inform the Manufacturer's Quality Service.

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