




## Endoscopic Biopsy Forceps

The Endoscopic Biopsy Forceps are available in both single-bite (without spike) and double-bite (with spike) design. The double-bite design offers the benefit of tissue impalement enabling multiple sampling in one procedure. The presence of a spike also provides a clear direction for lesion sampling and deeper biopsy penetration.

- Working length of 150 and 220cm
- Compatible with working channel of 2.8mm
- Sheath diameter of 2.4mm
- Single-bite (without spike) forceps for simple mono-biopsy sampling
- Double-bite (spike) forceps offer:
  - Multiple sampling
  - Tissue impalement
  - Directed lesion sampling
  - Deeper biopsy penetration



## Ordering information

For further information, or to place an order please call us on **01844 340 620**  
or email us at [info@aprmedtech.com](mailto:info@aprmedtech.com)

Product code	Description	Type	Working length (cm)	Sheath	Working channel
WS-2415BWH	Endoscopic Biopsy Forceps	With spike	150	2.4mm	2.8mm
WS-2415BTH	Endoscopic Biopsy Forceps	Without spike	150	2.4mm	2.8mm
WS-2422BWH	Endoscopic Biopsy Forceps	With spike	220	2.4mm	2.8mm
WS-2422BTH	Endoscopic Biopsy Forceps	Without spike	220	2.4mm	2.8mm

## Complementary products

We also offer a wide range of plastic stents for the management of occlusions within the biliary tree.

Description	Size (Fr)	Length (cm)	Box qty
PEBAX® Standard Biliary Stent	7, 8.5, 10, 11.5	5, 7, 10,12,15	10
Teflon Biliary Stent	7, 8.5, 10, 11.5	5, 7, 10,12,15	10
Amsterdam Biliary Stent	7, 8.5, 10, 11.5	5, 7, 10,12,15	10
Curved Pancreatic Stent	5, 7, 8.5, 11.5	3, 5, 7, 9, 10, 12, 15	10
Pigtail Pancreatic Stent	5, 7, 8.5, 11.5	3, 5, 7, 9, 10, 12, 15	10

[aprmedtech.com](http://aprmedtech.com)

@aprmedtech

APR Medtech Limited

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)

