


Accoat Guide Wires



About Accoat Guide Wires

HIGH QUALITY MANUFACTURE

Accoat Guide Wires are manufactured at the state-of-the-art production facilities of SP Medical, Denmark. They have more than 30 years' experience and expertise in the design and development of guide wires for the global medical device industry. If you would like to learn more a short video showing what goes in to making Accoat Guide Wires is available to view on our website.

ENGINEERED FOR YOU

Working with SP Medical allows us to offer you custom guide wire lengths and tip types, to suit specific clinical requirements. These can be manufactured in low volumes for individual clinicians or hospital departments. [Just let us know what length guide wire you require and we will do the rest.](#)

SERVICE AND SUPPORT

APR Medtech is a specialist independent medical technology company. We supply the Accoat range of guide wires to the NHS and private healthcare sector.

One of the key elements of our customer support service is product training. We want to be sure that our products are always used in the optimal way to the benefit of our customers and their patients.

If you would like one of our team to visit, either to attend a procedure or run a training session on one of our products, then please get in touch.

GET IN TOUCH

T. 01844 340 620

E. info@aprmedtech.com

Additional information is also available on our website: www.aprmedtech.com

Follow us on Twitter [@aprmedtech](https://twitter.com/aprmedtech)



Accoat Guide Wire Range


The Accoat range includes both [Super Stiff Seldinger](#) and [Lunderquist style](#) guide wires.

Each guide wire consists of a stainless steel core with PTFE coating and is available in various lengths with the option of a flexible straight or J-tip, making it well suited to a variety of interventional radiology and urology procedures.

RADIOPACITY: Stainless steel provides excellent radiopacity for clear visualisation under x-ray

CUSTOM LENGTHS: Available in 0.035" diameter and custom lengths to suit specific clinical requirements

PTFE COATING: Designed to minimise friction and aid tracking



ORDERING INFORMATION

Listed below are some of the popular sizes.

To enquire about custom length wires, please contact your local representative or a member of our customer support team.

Accoat Super Stiff Guide Wire


Product code	Description	Length (cm)	Diameter
HXX-S-C-035-80	Accoat Super Stiff Guide Wire Straight Tip 80cm	80	0.035"
HXX-S-C-035-150	Accoat Super Stiff Guide Wire Straight Tip 150cm	150	0.035"
HXX-S-C-035-260	Accoat Super Stiff Guide Wire Straight Tip 260cm	260	0.035"
HXX-J3-C-035-80	Accoat Super Stiff Guide Wire J-Tip 80cm	80	0.035"
HXX-J3-C-035-150	Accoat Super Stiff Guide Wire J-Tip 150cm	150	0.035"
HXX-J3-C-035-260	Accoat Super Stiff Guide Wire J-Tip 260cm	260	0.035"

Accoat Lunderquist Guide Wire*

Product code	Description	Length (cm)	Diameter
PTC-S-C-F10-035-080	Accoat Lunderquist 80cm Straight 10cm Flexible Tip	80	0.035"
PTC-S-C-F10-035-180	Accoat Lunderquist 180cm Straight 10cm Flexible Tip	180	0.035"
PTC-S-C-F10-035-260	Accoat Lunderquist 260cm Straight 10cm Flexible Tip	260	0.035"

* **DISTAL FLEXIBLE CM LENGTH:** Available in 5 (F5), 7 (F7), 10 (F10), 12 (F12), 15 (F15), and 17 (F17)

For additional information, or to order one of our products, please call us on [01844 340 620](tel:01844340620) or email info@aprmedtech.com



aprmedtech.com
[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX
t. 01844 340 620 f. 01844 340 621 e. info@aprmedtech.com

Manufactured by:

SP Medical