

Accoat™ Guide Wire

Consists of a stainless steel core with or without a PTFE coat. It is available with a flexible 4.5cm straight or J-tip. The J-tip radius is available in 1.5, 3, 6 or 15mm size. There are varying degrees of shaft stiffness available and wire length options vary from 30-450cm.



- Stainless steel provides excellent radiopacity for clear visualisation under x-ray
- PTFE coating designed to minimise friction and aid tracking
- Standard, Stiff, Extra Stiff and Super Stiff options offer an excellent platform choice for operator
- A 4.5cm flexible J-tip is provided with a radius of 1.5, 3, 6, or 15mm. A 4.5cm flexible straight tip option is also available




Accoat Guide Wire

Ordering information

For additional information, or to order this product, please call us on **01844 340 620** or email info@aprmedtech.com.

Listed below are some of the popular sizes.
Other sizes are available upon request.

| Product code | Description | Length (cm) | Diameter | Box qty |
|--------------------|---|-------------|----------|---------|
| J3-C-035-150 | Standard shaft, J-tip with 3mm radius, PTFE coated | 150 | 0.035" | 10 |
| J3-C-035-180 | Standard shaft, J-tip with 3mm radius, PTFE coated | 180 | 0.035" | 10 |
| J3-C-035-260 | Standard shaft, J-tip with 3mm radius, PTFE coated | 260 | 0.035" | 10 |
| H-X-X-J3-C-035-80 | Super stiff shaft, J-tip with 3mm radius, PTFE coated | 80 | 0.035" | 10 |
| H-X-X-J3-C-035-150 | Super stiff shaft, J-tip with 3mm radius, PTFE coated | 150 | 0.035" | 10 |
| H-X-X-J3-C-035-180 | Super stiff shaft, J-tip with 3mm radius, PTFE coated | 180 | 0.035" | 10 |
| H-X-X-J3-C-035-260 | Super stiff shaft, J-tip with 3mm radius, PTFE coated | 260 | 0.035" | 10 |



apr  medtech



aprmedtech.com

Manufactured by:

SP Medical

Additional product information is available at:

aprmedtech.com

APR Medtech Ltd

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX

t. +44 (0)1844 340 620

f. +44 (0)1844 340 621

e. info@aprmedtech.com