


# Poseidon™ Hydrophilic Nitinol Guide Wire

One of the key outstanding features of the Poseidon Hydrophilic Nitinol Guide Wire is its long lasting hydrophilic coating, which has been designed to minimise friction and aid tracking. The Poseidon is also available in a variety of tip types, lengths and degrees of shaft stiffness. The wide range provides options to support most procedural eventualities.

- Excellent radiopacity for clear visualisation under x-ray
- Lubricious long-lasting hydrophilic polyurethane coating designed to minimise friction and aid tracking
- Nitinol core offering kink resistance and excellent flexibility
- Various tips types and shaft stiffness options available



## Ordering information

For further information or to place an order, please call us on **01844 340 620** or email us at [info@aprmtech.com](mailto:info@aprmtech.com).


Listed below are some of the popular sizes. Other sizes are available upon request.

Product code	Description	Length (cm)	Diameter	Box qty
HY-N-A-ST-035-D-150	Hydrophilic, angled tip, standard tip stiffness, standard shaft	150	0.035"	5
HY-N-A-ST-035-D-180	Hydrophilic, angled tip, standard tip stiffness, standard shaft	180	0.035"	5
HY-N-A-ST-035-D-260	Hydrophilic, angled tip, standard tip stiffness, standard shaft	260	0.035"	5
HY-N-A-ST-035-H-150	Hydrophilic, angled tip, standard tip stiffness, stiff shaft	150	0.035"	5
HY-N-A-ST-035-H-180	Hydrophilic, angled tip, standard tip stiffness, stiff shaft	180	0.035"	5
HY-N-A-ST-035-H-260	Hydrophilic, angled tip, standard tip stiffness, stiff shaft	260	0.035"	5

## Product code identifier

HY— N — S/A/R — ST/SO/F — 0xx — D/H — L  
 (example: HY-N-A-ST-035-D-150)

HY	Hydrophilic coating
N	Nitinol core
S	Straight tip
A	Angled tip
R	Re-shaping tip
ST	Standard tip stiffness
SO	Soft tip stiffness
F	Floppy tip stiffness
0xx	Diameter of the guidewire 0.018" 0.025" 0.032" 0.035" 0.038"
H	Stiff shaft
D	Standard shaft
L	Length of the guidewire (from 30cm—260cm)



[aprmtech.com](http://aprmtech.com)  
[@aprmtech](https://twitter.com/aprmtech)

APR Medtech Limited

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
 t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmtech.com](mailto:info@aprmtech.com)

