

The Salgado™ Driver

For accurate and secure placement of the
Haemorrhage Occluder Pin



For use in
combination
with the
Haemorrhage
Occluder Pin

A quick and simple to use device whereby a gentle tap of the driver handle delivers adequate force to fully insert the Haemorrhage Occluder Pin into the sacrum.

- Reusable
- Curved section rotates around the pubic symphysis
- Minimal force required to quickly drive the HOP into the sacrum

The Salgado™ Driver

For accurate and secure placement of the Haemorrhage Occluder Pin


Ordering information

For additional information, or to order this product, please call our Customer Support Team on [01844 340 620](tel:01844340620) or email info@aparmedtech.com.

Code	Description	Box qty
CR1007	Haemorrhage Occluder Pin Standard - 10mm	2
CR1014	Haemorrhage Occluder Pin Large - 14mm	2
CR2000	Salgado Driver Standard - 23cm length	1
CR2005	Salgado Driver Large - 35cm length	1

Did you know...

The Haemorrhage Occluder Pin is available in 10mm & 14mm to tackle bleeds of varying size and location.



Manufactured by:



© 2016 APR Medtech Ltd. Hemorrhage Occluder Pin and Salgado Driver are trademarks of Surgin, Inc. All rights reserved.

SUR-SALGADO-A5-v1-0716. Date of preparation: July 2016.

Additional product information is available at:

aparmedtech.com

apr  medtech



aparmedtech.com

APR Medtech Ltd

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX

t. +44 (0)1844 340 620

f. +44 (0)1844 340 621

e. info@aparmedtech.com