



# Yangtze<sup>®</sup> angioplasty balloon catheters

The Yangtze<sup>®</sup> range of angioplasty balloon catheters are available in a broad range of sizes and diameters ideal for a wide range of coronary interventions.

## Yangtze<sup>®</sup> $\mu$


semi-compliant balloon catheter



- Low profile (0.016") to aid advancement through narrowed, heavily calcified lesions
- Length: 10 to 30 mm | Diameter: 1.5 to 5.00 mm
- Nominal pressure: 6 atm | Rated Burst Pressure: 16 atm
- Guidewire compatibility: 0.014"
- Guiding catheter compatibility: 5F for Ø1.50 to 4.00 mm, 6F for Ø4.50 and 5.00 mm

## Yangtze<sup>®</sup> NC

non-compliant balloon catheter



- Length: 10, 15 and 20 mm | Diameter: 2.00 to 4.00 mm
- Nominal Pressure (NP): 12 atm
- Rated Burst Pressure (RBP): 20 atm
- Average Burst Pressure (ABP): 27 atm
- Guidewire compatibility: 0.014"
- Guiding catheter compatibility: 5F

## Yangtze<sup>®</sup> $\mu^{sb}$

bifurcation side branch balloon catheter



- Length: 10mm | Diameter: 2.00 and 3.00 mm
- Nominal pressure: 6 atm
- Rated Burst Pressure: 16 atm
- Guidewire compatibility: 0.014"
- Guiding catheter compatibility: 5F

## Ordering information

For further information or to place an order, please call us on 01844 340 620 or email us at [info@aprmtech.com](mailto:info@aprmtech.com).

### Yangtze® $\mu$ semi-compliant balloon catheter

Length (mm)	Diameter (mm)						
	1.00	1.10	1.25	1.50	2.00	2.25	2.50
8	PTMU0054	PTMU0055	PTMU0056	-	-	-	-
10	PTMU0057	PTMU0058	PTMU0059	PTMU0001	PTMU0002	PTMU0003	PTMU0004
15	PTMU0060	PTMU0061	PTMU0062	PTMU0010	PTMU0011	PTMU0012	PTMU0013
20	PTMU0063	PTMU0064	PTMU0065	PTMU0019	PTMU0020	PTMU0021	PTMU0022
25	PTMU0066	PTMU0067	PTMU0068	-	PTMU0028	PTMU0029	PTMU0030
30	-	-	-	-	PTMU0036	PTMU0037	PTMU0038

Length (mm)	Diameter (mm)						
	2.75	3.00	3.25	3.50	4.00	4.50	5.00
8	-	-	-	-	-	-	-
10	PTMU0005	PTMU0006	PTMU0007	PTMU0008	PTMU0009	PTMU0044	PTMU0045
15	PTMU0014	PTMU0015	PTMU0016	PTMU0017	PTMU0018	PTMU0046	PTMU0047
20	PTMU0023	PTMU0024	PTMU0025	PTMU0026	PTMU0027	PTMU0048	PTMU0049
25	PTMU0031	PTMU0032	PTMU0033	PTMU0034	PTMU0035	PTMU0050	PTMU0051
30	PTMU0039	PTMU0040	PTMU0041	PTMU0042	PTMU0043	PTMU0052	PTMU0053

### Yangtze® NC non-compliant balloon catheter

Length (mm)	Diameter (mm)				
	2.0	2.5	3.0	3.5	4.0
10	PTCP0001	PTCP0002	PTCP0003	PTCP0004	PTCP0005
15	PTCP0007	PTCP0008	PTCP0009	PTCP0010	PTCP0011
20	PTCP0013	PTCP0014	PTCP0015	PTCP0016	PTCP0017

### Yangtze® $\mu^{sb}$ bifurcation side branch balloon catheter

Length (mm)	Diameter (mm)		
	2.0	2.5	3.0
10	PTMU0001	PTMU0002	PTMU0003

[aprmtech.com](http://aprmtech.com)

@aprmtech

APR Medtech Limited

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX

t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmtech.com](mailto:info@aprmtech.com)

