




Devoir

Sirolimus Coated Balloon Catheter

- Bio-sphere drug delivery technology
- Low rate of MACE with safety and efficacy reported up to 12 months*
- 5Fr compatible for all sizes
- Drug between balloon folds for reduced in-transit loss of drug (10%)

Clinical data*

12 months clinical results (277 patients)



*Data on file – clinical data all-comers population.



Devoir

Sirolimus Coated Balloon Catheter

Product specification

Drug	Sirolimus
Drug dose	1.27 µg / mm ²
Drug carrier	Phospholipid based excipient
Balloon material	Polyamide blend (semi-compliant)
Catheter design	Rapid exchange
Shaft diameter - prox.	1.7F
Shaft diameter - dist.	2.5F
Usable catheter length	140cm
Tip profile	0.016"
Nominal Pressure (NP)	6 bar
Rated Burst Pressure (RBP)	16 bar (14 bar for 4.00 / 25 to 40mm)
Guide cath. compatibility	5Fr (0.056" ID) for all sizes
Guidewire compatibility	0.014" maximum recommended

Ordering information

Diameter (mm)	Length						
	10	15	20	25	30	35	40
1.50	CDE15010	CDE15015	CDE15020	CDE15025	CDE15030	CDE15035	CDE15040
2.00	CDE20010	CDE20015	CDE20020	CDE20025	CDE20030	CDE20035	CDE20040
2.25	CDE22510	CDE22515	CDE22520	CDE22525	CDE22530	CDE22535	CDE22540
2.50	CDE25010	CDE25015	CDE25020	CDE25025	CDE25030	CDE25035	CDE25040
2.75	CDE27510	CDE27515	CDE27520	CDE27525	CDE27530	CDE27535	CDE27540
3.00	CDE30010	CDE30015	CDE30020	CDE30025	CDE30030	CDE30035	CDE30040
3.50	CDE35010	CDE35015	CDE35020	CDE35025	CDE35030	CDE35035	CDE35040
4.00	CDE40010	CDE40015	CDE40020	CDE40025	CDE40030	CDE40035	CDE40040

For additional information, or to order this product, please call us on **01844 340 620** or email info@aprmedtech.com

Additional product information is available at:
aprmedtech.com

apr  medtech



aprmedtech.com

In partnership with:



Manufactured by Envision Scientific Pvt. Ltd.

APR Medtech Ltd

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX

t. +44 (0)1844 340 620

f. +44 (0)1844 340 621

e. info@aprmedtech.com