


Interventional Radiology



OUR COMPANY

APR Medtech is a specialist independent medical technology company providing high quality medical devices and support services to the NHS and private healthcare sector.

Our products offer intelligent clinical solutions designed to help advance healthcare and improve patient outcomes. Our extensive experience within the healthcare industry, deep knowledge of the clinical areas in which we work and a genuine passion for what we do, ensures the delivery of real value to our customers.

SERVICE AND SUPPORT

One of the key elements of our customer support service is product training. We want to be sure that our products are always used in the optimal way to the benefit of both our customers and their patients.


If you would like one of our team to visit, either to attend a procedure or run a training session on one of our products, then please get in touch. We will be more than happy to oblige.

GET IN TOUCH

For additional information please call us on [01844 340 620](tel:01844340620) or email info@aprmedtech.com.

Additional information is also available on our website: www.aprmedtech.com.

Follow us on Twitter [@aprmedtech](https://twitter.com/aprmedtech) 



MULTIPURPOSE DRAINAGE

MPD - Multipurpose Drainage Catheter Sets

Our MPD Multipurpose Drainage Catheters and Sets offer a new solution to the challenge of maintaining a soft, kink free material with the need to prevent 'ruckle-back' during insertion.

SOFT, KINK FREE catheter material improves patient comfort and prevents unnecessary catheter replacement caused by harder, 'kink resistant' materials.

SECURE LOCKING SYSTEM to give reliable pigtail locking without cumbersome 'keys' or bulky hubs. Self-adhesive removal instructions (for addition to the patient's notes) are provided with every locking drain.

EASYTIP™ CORE to prevent 'ruckle-back', minimising resistance during insertion; reducing discomfort and speeding correct placement (8.4F-14F).

COMPLETE WITH DRAINAGE CONNECTION SET to reduce costs and make our MPD Catheters even better value for your department or ward.


ORDERING INFORMATION

Additional French sizes and lengths are available upon request.

Description	Product code	Size	Length (CM)	Box qty
Multipurpose Drainage Catheter - Non Locking Contains pigtail catheter, trocar and stiffener set and drainage bag connection set with 30cm tubing	R58810-06-30	6F	30	10
	R58810-84-30	8F	30	10
	R58810-10-30	10F	30	10
	R58810-14-30	14F	30	10
Multipurpose Drainage Catheter - Locking Contains pigtail catheter, trocar and stiffener set and drainage bag connection set with 30cm tubing	R58830-06-30	6F	30	10
	R58830-84-30	8F	30	10
	R58830-10-30	10F	30	10
	R58830-14-30	14F	30	10
MPD Exchange Pigtail Catheter	R58851-84-30	8F	30	10
Drainage Connection Set Funnel end connector, 30cm tubing and luer stopcock with rotating locking ring	R58915	N/A	30	10

Sterile. Single patient use only.

Includes
Drainage
Connection
Set



MULTIPURPOSE DRAINAGE PACKS

Ascitex™ Multipurpose Drainage Procedure Pack

The Ascitex™ Multipurpose Drainage Procedure Pack contains all the items necessary to safely and conveniently complete the procedure in the department or ward. It is indicated for percutaneous introduction of a drain for paracentesis, the drainage of pleural fluid or similar cavity drainage.

MULTIPURPOSE DRAIN with soft, kink resistant catheter material, improves patient comfort and prevents unnecessary catheter replacement caused by harder materials.

NITINOL GUIDEWIRE with soft, non-traumatic platinum tip to ensure safe delivery of the catheter into the abdominal cavity or organ.

SECURE LOCKING SYSTEM to give reliable pigtail locking without cumbersome 'keys' or bulky hubs. Self-adhesive removal instructions (for addition to the patient's notes) are provided with every locking drain.

SELF-ADHESIVE DRAIN FIXATION is included in every kit to ensure comfort and extra security.

COMPLETE PROCEDURE TRAY including syringes, needles, drapes, swabs and ascites drainage bag with tap.



Includes
safety scalpel

ORDERING INFORMATION

Description	Product code	Box qty
Ascitex™ Multipurpose Drainage Procedure Pack	R58800-84-30	Supplied in single units

Sterile. Single patient use only.

Pack contains:

- 8.4F locking pigtail catheter
- 6 holes x 30cm, 20ml
- 12ml syringes, luer lock
- 21G x 1.5" needle x 2
- 25G x 5/8" needle
- Scalpel No.11
- 18G x 7cm Seldinger needle
- Nitinol guidewire 0.035" x 80cm x 6cm soft tip
- 9F x 8in dilator
- 30cm drain bag connection set + 2 way tap
- 2000ml drainage bag with drain port
- 60ml gallipot, 10cm² x-ray detectable gauze swabs x 5
- Drain Guard fixation device
- 50 x 60m adhesive fenestrated drape
- 100 x 100cm sterile field wrap
- Locking pigtail insertion and removal instructions

INTRODUCER NEEDLES

Introducer Needles


Two and three-part introducer needle sets for single puncture percutaneous access to allow placement of up to a 0.038" OD wire guide.

LOCKING HUB all introducer needles feature a fully locking hub system to make insertion easier and more secure.


RADIOPAQUE & CLEAR 'Kellet style' versions are available as standard to suit all users. Outer permits up to 0.038" wire, inner cannula up to 0.032" wire.

SMOOTH TRANSITION between the tri-facet or Chiba needle, cannula and outer sheath to ensure trouble-free insertion with the minimum of trauma.

SILICONISED needle sets ensure smooth, drag free withdrawal after insertion.



Fully locking hubs



DUO™ Chiba bevel

ORDERING INFORMATION

Description	Product code	Sheath	Length (CM)	Box qty
Introducer Needle Three part, 5F outer sheath, 19G inner cannula tube, 21G tri-facet trocar needle.	R58900	R/O	19.5	10
	R58900-00-KL	Clear	19.5	10
	R58900-00-15	R/O	15	10
	R58900-KL-15	Clear	15	10
DUO™ Introducer Needle Two part, 5F outer sheath, 19G inner Chiba needle	R58900-CH-15	Clear	15	10

Sterile. Single patient use only.

RADIOPAQUE DILATORS & ACCESS SET

Radiopaque Dilators


FULL SIZE RANGE 5F to 16F with colour coding for easy identification.

LUER HUB for secure connection.

SMOOTH TRANSITION from tip to guidewire to ensure easy access.

Description	French size	Product Code	Box qty
Radiopaque Dilators Luer lock, polyethylene, barium-loaded, tapered dilators All sizes are 21cm in length	5	R58921-05-21	10
	6	R58921-06-21	
	7	R58921-07-21	
	8	R58921-08-21	
	9	R59821-09-21	
	10	R58921-10-21	
	12	R58921-12-21	
	14	R58921-14-21	
	16	R58921-16-21	

Sterile. Single patient use only



Percutaneous Access Set

A convenient kit for percutaneous access which includes:

0.018" PTFE J-tip guidewire.


0.035" PTFE J-tip guidewire.

4, 5, 6, 7 and 8F radiolucent dilators.

18G introducer needle.

Description	Product Code	Box qty
Percutaneous Access Set 0.018" PTFE J-tip guidewire 80cm 0.035" PTFE J-tip guidewire 80cm 4,5,6,7 and 8F radiolucent dilators, 21cm length 18G Introducer needle x 7cm	R58912	10

Sterile. Single patient use only.




DRAINAGE ACCESSORIES

Drainage Connection Set

Includes funnel end connector, luer stopcock with rotating locking ring and 30cm tubing.

Description	Product code	Box qty
Drainage Connection Set	R58915	10

Note: this connection set is included with each MPD Catheter Set.



For additional information, or to order one of our products, please call us on 01844 340 620 or email info@aprimedtech.com.

Further product information and support materials are also available on our website at www.aprimedtech.com.



For additional information, or to order one of our products, please call us on 01844 340 620 or email info@aprmedtech.com.

aprmedtech.com
[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX
t. +44 1844 340 620 f. +44 1844 340 621 e. info@aprmedtech.com



Rocket Medical PLC
Rocket Medical plc., Sedling Road,
Washington Tyne & Wear, NE38 9BZ



Rocketmedical