


Ascitex™ Procedure Pack
for day-case paracentesis



OUR COMPANY

APR Medtech is a specialist independent medical technology company providing high quality medical devices and support services to the NHS and private healthcare sector.

Our products offer intelligent clinical solutions designed to help advance healthcare and improve patient outcomes. Our extensive experience within the healthcare industry, deep knowledge of the clinical areas in which we work and a genuine passion for what we do, ensures the delivery of real value to our customers.

SERVICE AND SUPPORT

One of the key elements of our customer support service is product training. We want to be sure that our products are always used in the optimal way to the benefit of both our customers and their patients.

If you would like one of our team to visit, either to attend a procedure or run a training session on one of our products, then please get in touch. We will be more than happy to oblige.

GET IN TOUCH

For additional information please call us on [01844 340 620](tel:01844340620) or email info@aprmedtech.com

Additional information is also available on our website: www.aprmedtech.com

Follow us on Twitter [@aprmedtech](https://twitter.com/aprmedtech) 



PARACENTESIS

Ascitex™ Multipurpose Drainage Procedure Pack

The Ascitex™ Multipurpose Drainage Procedure Pack contains all the items necessary to safely and conveniently complete the procedure in the department or ward. It is indicated for percutaneous introduction of a drain for paracentesis, the drainage of pleural fluid or similar cavity drainage.

MULTIPURPOSE DRAIN with soft, kink resistant catheter material, improves patient comfort and prevents unnecessary catheter replacement caused by harder materials.

NITINOL GUIDEWIRE with soft, non-traumatic platinum tip to ensure safe delivery of the catheter into the abdominal cavity or organ.

SECURE LOCKING SYSTEM to give reliable pigtail locking without cumbersome 'keys' or bulky hubs. Self-adhesive removal instructions (for addition to the patient's notes) are provided with every locking drain.

SELF-ADHESIVE DRAIN FIXATION is included in every kit to ensure comfort and extra security.

COMPLETE PROCEDURE TRAY including syringes, needles, drapes, swabs and ascites drainage bag with tap.



Includes
safety scalpel

ORDERING INFORMATION

Description	Product code	Box qty
Ascitex™ Multipurpose Drainage Procedure Pack	R58800-84-30	Supplied in single units

Sterile. Single patient use only.

Pack contents:

- 10ml luer lock syringe
- 20ml luer lock syringe
- 25G x 16mm (5/8") needle
- 21G x 38mm (1.5") needle
- 18G x 70mm Seldinger needle
- Safety scalpel (No.11)
- Nitinol guidewire 0.035" x 80cm with 6cm soft tip
- 9F x 21cm (8") dilator
- 8.4F x 30cm locking pigtail catheter
- Drain Guard fixation device
- 30cm drainage bag connection set with 2-way tap
- 2000ml drainage bag with drain port
- 60ml gallipot
- 5x 10cm² X-ray detectable gauze swabs
- 50 x 60m adhesive fenestrated drape
- 100 x 100cm sterile field wrap
- Locking pigtail catheter insertion instructions
- Locking pigtail catheter removal instructions sticker

For additional information, or to order one of our products, please call us on
01844 340 620 or email info@aprmedtech.com

aprmedtech.com
[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX
t. +44 1844 340 620 f. +44 1844 340 621 e. info@aprmedtech.com



Rocket Medical PLC
Rocket Medical plc., Sedling Road,
Washington Tyne & Wear, NE38 9BZ



Rocketmedical