

Manawa
Ureteral Access Sheath



OUR COMPANY

APR Medtech is a specialist independent medical technology company providing high quality medical devices and support services to the NHS and private healthcare sector.

Our products offer intelligent clinical solutions designed to help advance healthcare and improve patient outcomes. Our extensive experience within the healthcare industry, deep knowledge of the clinical areas in which we work and a genuine passion for what we do, ensures the delivery of real value to our customers.

SERVICE AND SUPPORT

One of the key elements of our customer support service is product training. We want to be sure that our products are always used in the optimal way to the benefit of both our customers and their patients.

If you would like one of our team to visit, either to attend a procedure or run a training session on one of our products, then please get in touch. We will be more than happy to oblige.

GET IN TOUCH

For additional information please call us on [01844 340 620](tel:01844340620) or email info@aprmedtech.com

Additional information is also available on our website: www.aprmedtech.com

Follow us on Twitter [@aprmedtech](https://twitter.com/aprmedtech) 




Ureteral Access Sheath

Provides ureteral dilation with a continuous working channel for the introduction of ureteroscopes and instruments during ureteroscopic procedures.

1—KINK-RESISTANT SHEATH: made with radiopaque stainless steel coil for improved stability and trackability.

2—RADIOPAQUE RING MARKER: facilitates placement of the sheath tip under x-ray.

3—FLEXIBLE TAPERED DILATOR TIP: ensures gentle atraumatic insertion and dilatation.



4—ERGONOMIC HUB HANDLE: for confident and easy manipulation.

5—PROPRIETARY HYDROPHILIC COATING: applied to the sheath and dilator tip to provide flawless insertion.

6—SEAMLESS TRANSITION: facilitates insertion of the sheath and reduces risk of trauma

Ureteral Access Sheath

ORDERING INFORMATION

Description	Product code	Internal diameter French Size	External diameter French size	Length
Manawa Access Sheath	UAS-0109520S	9.5	11.5	20cm
Manawa Access Sheath	UAS-0109528S	9.5	11.5	28cm
Manawa Access Sheath	UAS-0109535S	9.5	11.5	35cm
Manawa Access Sheath	UAS-0109545S	9.5	11.5	45cm
Manawa Access Sheath	UAS-0110720S	10.7	12.7	20cm
Manawa Access Sheath	UAS-0110728S	10.7	12.7	28cm
Manawa Access Sheath	UAS-0110735S	10.7	12.7	35cm
Manawa Access Sheath	UAS-0110745S	10.7	12.7	45cm
Manawa Access Sheath	UAS-0112020S	12.0	14.0	20cm
Manawa Access Sheath	UAS-0112028S	12.0	14.0	28cm
Manawa Access Sheath	UAS-0112035S	12.0	14.0	35cm
Manawa Access Sheath	UAS-0112045S	12.0	14.0	45cm
Manawa Access Sheath	UAS-0114020S	14.0	16.0	20cm
Manawa Access Sheath	UAS-0114028S	14.0	16.0	28cm
Manawa Access Sheath	UAS-0114035S	14.0	16.0	35cm
Manawa Access Sheath	UAS-0114045S	14.0	16.0	45cm

For additional information, or to order one of our products, please call us on 01844 340 620 or email info@aprmedtech.com

aprmedtech.com
[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX
t. +44 1844 340 620 f. +44 1844 340 621 e. info@aprmedtech.com



Natec Medical Limited
Maeva Centre Building,
Silicon Avenue, Ebene Business Park,
Réduit 72 201 Mauritius

