


Manawa  
Ureteral Access Sheath



# URETERAL ACCESS SHEATH


## Manawa Ureteral Access Sheath

Provides ureteral dilation with a continuous working channel for the introduction of ureteroscopes and instruments during ureteroscopic procedures.

**1—KINK-RESISTANT SHEATH:** made with radiopaque stainless steel coil for improved stability and trackability

**2—RADIOPAQUE RING MARKER:** facilitates placement of the sheath tip under x-ray

**3—FLEXIBLE TAPERED DILATOR TIP:** ensures gentle atraumatic insertion and dilatation



**4—SEAMLESS TRANSITION:** facilitates insertion of the sheath and reduces risk of trauma

**5—PROPRIETARY HYDROPHILIC COATING:** applied to the sheath and dilator tip to provide flawless insertion

**6—ERGONOMIC HUB HANDLE:** for confident and easy manipulation

# URETERAL ACCESS SHEATH

## Manawa Ureteral Access Sheath

### ORDERING INFORMATION

For additional information, or to order one of our products, please call us on [01844 340 620](tel:01844340620) or email [info@aprmedtech.com](mailto:info@aprmedtech.com)

Description	Product code	Internal diameter French Size	External diameter French size	Length
Manawa Access Sheath	UAS-0109520S	9.5	11.5	20cm
Manawa Access Sheath	UAS-0109528S	9.5	11.5	28cm
Manawa Access Sheath	UAS-0109535S	9.5	11.5	35cm
Manawa Access Sheath	UAS-0109545S	9.5	11.5	45cm
Manawa Access Sheath	UAS-0110720S	10.7	12.7	20cm
Manawa Access Sheath	UAS-0110728S	10.7	12.7	28cm
Manawa Access Sheath	UAS-0110735S	10.7	12.7	35cm
Manawa Access Sheath	UAS-0110745S	10.7	12.7	45cm
Manawa Access Sheath	UAS-0112020S	12.0	14.0	20cm
Manawa Access Sheath	UAS-0112028S	12.0	14.0	28cm
Manawa Access Sheath	UAS-0112035S	12.0	14.0	35cm
Manawa Access Sheath	UAS-0112045S	12.0	14.0	45cm
Manawa Access Sheath	UAS-0114020S	14.0	16.0	20cm
Manawa Access Sheath	UAS-0114028S	14.0	16.0	28cm
Manawa Access Sheath	UAS-0114035S	14.0	16.0	35cm
Manawa Access Sheath	UAS-0114045S	14.0	16.0	45cm

For additional information, or to order one of our products, please call us on  
01844 340 620 or email [info@aprmedtech.com](mailto:info@aprmedtech.com)

[aprmedtech.com](http://aprmedtech.com)  
[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited  
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)



Natec Medical Limited  
Maeva Centre Building,  
Silicon Avenue, Ebene Business Park,  
Réduit 72 201 Mauritius

