


## Tahina

### Nephrostomy Balloon Catheter



## OUR COMPANY

APR Medtech is a specialist independent medical technology company providing high quality medical devices and support services to the NHS and private healthcare sector.

Our products offer intelligent clinical solutions designed to help advance healthcare and improve patient outcomes. Our extensive experience within the healthcare industry, deep knowledge of the clinical areas in which we work and a genuine passion for what we do, ensures the delivery of real value to our customers.

## SERVICE AND SUPPORT

One of the key elements of our customer support service is product training. We want to be sure that our products are always used in the optimal way to the benefit of both our customers and their patients.

If you would like one of our team to visit, either to attend a procedure or run a training session on one of our products, then please get in touch. We will be more than happy to oblige.

## GET IN TOUCH

For additional information please call us on [01844 340 620](tel:01844340620) or email [info@aprmedtech.com](mailto:info@aprmedtech.com)

Additional information is also available on our website: [www.aprmedtech.com](http://www.aprmedtech.com)

Follow us on Twitter [@aprmedtech](https://twitter.com/aprmedtech) 



# NEPHROSTOMY BALLOON CATHETER

## Tahina Nephrostomy Balloon Catheter

Powerful and versatile balloon dilator making establishment of a nephrostomy tract a safer and faster procedure.

**ULTRA LOW BALLOON PROFILE** to enhance access for difficult stone impacted calyces


**PROPRIETARY HYDROPHILIC COATING** applied to the balloon surface and tip eases catheter advancement

**SHORT TAPERED TIP DESIGN** and short tapered balloon shoulders to enable close dilatation near the calculi

**HIGH RADIAL DILATATION BALLOON** engineered to withstand high pressure for confident dilatation of the nephrostomy tract and convenient placement of the sheath

**HIGHLY RADIOPAQUE** marker band and PTFE sheath to facilitate accurate balloon placement

Includes  
radiopaque  
PTFE sheath with  
bevelled end



### ORDERING INFORMATION

Description	Product code	French Size	Guide wire size	Catheter length	Inflated diameter	Balloon length	Inflation pressure	Box qty
Tahina Nephrostomy Balloon Catheter 24Fr	UDB-0308120S	6.8F	0.038"	55cm	24Fr	120cm	17atm	1
Tahina Nephrostomy Balloon Catheter 30Fr	UDB-0310120S	6.8F	0.038"	55cm	30Fr	120cm	17atm	1

For additional information, or to order one of our products, please call us on  
01844 340 620 or email [info@aprmedtech.com](mailto:info@aprmedtech.com)

[aprmedtech.com](http://aprmedtech.com)  
[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited  
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)



Natec Medical Limited  
Maeva Centre Building,  
Silicon Avenue, Ebene Business Park,  
Réduit 72 201 Mauritius

