


NEXENT

BILIARY SELF-EXPANDING NITINOL STENT



Functional multilayer
coating technology

Novel radiopaque
marker




Excellent shape
conformity


Unique dispersed
hook and cross
structure




Excellent visibility



Shape conformity



Lumen conformity



Simple deployment

NEXENT

BILIARY STENT

Unique dispersed hook and cross structure

- Optimal balance of high radial and low axial force
- Excellent lumen conformity and compression resistance

Novel radiopaque marker design

- 12 flat disc-shaped tantalum markers
- Clear fluoroscopic visualisation


Protective silicone bilayer

- Encapsulates and protects stent markers
- Prevents tissue ingrowth and facilitates atraumatic removal

Highly uniform membrane thickness


- Reduces deployment force for smooth stent delivery

Available uncovered, partially covered or fully covered



EASY-TO-USE PULL BACK DELIVERY SYSTEM

- 8Fr TTS endoscopic (1800mm) and percutaneous (600mm) delivery system available
- Resheath up to three times



Manufactured by:



ORDERING INFORMATION

For additional information, please call us on 01844 340 620 or email info@aprmedtech.com

PRODUCT CODE			STENT		DELIVERY SYSTEM	
Uncovered	Partially Covered	Covered	Diameter (mm)	Length (mm)	Diameter (Fr)	Length (mm)
BNC-10-040-180-8	BPC-10-040-180-8	BCC-10-040-180-8	10	40	8	1800*
BNC-10-050-180-8	BPC-10-050-180-8	BCC-10-050-180-8		50		
BNC-10-060-180-8	BPC-10-060-180-8	BCC-10-060-180-8		70		
BNC-10-070-180-8	BPC-10-070-180-8	BCC-10-070-180-8		70		
BNC-10-080-180-8	BPC-10-080-180-8	BCC-10-080-180-8		80		
BNC-10-090-180-8	BPC-10-090-180-8	BCC-10-090-180-8		90		
BNC-10-100-180-8	BPC-10-100-180-8	BCC-10-100-180-8		100		
BNC-10-110-180-8	BPC-10-110-180-8	BCC-10-110-180-8		110		
BNC-10-120-180-8	BPC-10-120-180-8	BCC-10-120-180-8		120		

*Delivery system also available in 600mm length

Manufactured by:

