


NEXENT

COLORECTAL SELF-EXPANDING NITINOL STENT



Unique dispersed hook
and cross structure

Functional multilayer
coating technology




Novel radiopaque
marker


Excellent shape
conformity




Excellent visibility



Shape conformity



Lumen conformity



Simple deployment

NEXENT

COLORECTAL STENT

Unique dispersed hook and cross structure

- Optimal balance of high radial and low axial force
- Excellent lumen conformity and compression resistance

Novel radiopaque marker design

- 12 flat disc-shaped tantalum markers
- Clear fluoroscopic visualisation


Protective silicone bilayer

- Encapsulates and protects stent markers
- Prevents tissue ingrowth and facilitates atraumatic removal

Highly uniform membrane thickness


- Reduces deployment force for smooth stent delivery

Available uncovered, partially covered or fully covered



EASY-TO-USE PULL BACK DELIVERY SYSTEM

- 10.2Fr TTS delivery system
- Resheath up to three times



Manufactured by:



ORDERING INFORMATION

For additional information, please call us on 01844 340 620 or email info@aprmedtech.com

PRODUCT CODE			STENT		DELIVERY SYSTEM	
Uncovered	Partially Covered	Covered	Diameter (mm)	Length (mm)	Diameter (Fr)	Length (mm)
CNC-24-060-230-10	CPC-24-060-230-10	CCC-24-060-230-10	24	60	10.2	2300 *
CNC-24-070-230-10	CPC-24-070-230-10	CCC-24-070-230-10		70		
CNC-24-080-230-10	CPC-24-080-230-10	CCC-24-080-230-10		80		
CNC-24-090-230-10	CPC-24-090-230-10	CCC-24-090-230-10		90		
CNC-24-100-230-10	CPC-24-100-230-10	CCC-24-100-230-10		100		
CNC-24-110-230-10	CPC-24-110-230-10	CCC-24-110-230-10		110		
CNC-24-120-230-10	CPC-24-120-230-10	CCC-24-120-230-10		120		
CNC-24-130-230-10	CPC-24-130-230-10	CCC-24-130-230-10		130		
CNC-24-140-230-10	CPC-24-140-230-10	CCC-24-140-230-10		140		

*Delivery system also available in 1400mm length

Manufactured by:

