



# GATE BILIARY

## SELF-EXPANDING NITINOL STENT

**Woven stent construction** with small, unfixed cells offers high flexibility with minimal shortening

**The nitinol weave** rotates in both planes and in different directions for optimal conformability





**Gold radiopaque markers** provide excellent visualisation under fluoroscopy

 Simple to use **8F TTS delivery system**

 Available in **starting diameter of 8mm**




### Uncovered stent


Indicated for the management of malignant biliary strictures

-  **High flexibility** with minimal shortening
-  **Excellent conformability and flexibility** to accommodate tortuous anatomy

### Covered stent

Indicated for the management of malignant and benign biliary strictures

-  **Silicone coating** designed to resist tissue ingrowth and facilitate safe removal
-  **Retrieval string on distal end** to aid safe stent removal
-  Option of **flared or straight ends**



## DESIGN CHOICE

The Gate Biliary Stent is available in a broad range of sizes, with either an endoscopic or percutaneous delivery system. Please contact us with your preferred stent specification and we will provide you with additional ordering information.

### Endoscopic application

STENT TYPE	DELIVERY SYSTEM		STENT	
	Diameter (Fr)	Length (mm)	Diameter (mm)	Length (cm)
UNCOVERED / STRAIGHT	8	1800	8, 10, 12	4, 5, 6, 7, 8, 9, 10, 11, 12
COVERED / STRAIGHT	8	1800	8, 10	4, 5, 6, 7, 8, 9, 10, 11, 12
COVERED / FLARED	8	1800	8, 10	4, 5, 6, 7, 8, 9, 10, 11, 12

### Percutaneous application

STENT TYPE	DELIVERY SYSTEM		STENT	
	Diameter (Fr)	Length (mm)	Diameter (mm)	Length (cm)
UNCOVERED / STRAIGHT	8	500	8, 10, 12	4, 5, 6, 7, 8, 9, 10, 11, 12
COVERED / STRAIGHT	8	500	8, 10	4, 5, 6, 7, 8, 9, 10, 11, 12
COVERED / FLARED	8	500	8, 10	4, 5, 6, 7, 8, 9, 10, 11, 12

For additional information, please call us on [01844 340 620](tel:01844340620) or email [info@aprmedtech.com](mailto:info@aprmedtech.com)

Manufactured by:



[www.mediwood.co.kr](http://www.mediwood.co.kr)

Supplied exclusively in the UK by:



[aprmedtech.com](http://aprmedtech.com)  
@aprmedtech

APR Medtech Limited  
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)

