



# GATE COLORECTAL

## SELF-EXPANDING NITINOL STENT

**Woven stent construction**  
with small unfixed cells offers high flexibility with minimal shortening

**The nitinol weave** rotates in both planes and in different directions for optimal conformability

**Gold radiopaque markers** provide excellent visualisation under fluoroscopy



- Option of **flared or straight ends**
- Simple to use **10F TTS delivery system**
- Available with a **starting diameter of 18mm**

### Uncovered stent


Indicated for the management of malignant colorectal strictures



- High flexibility** with minimal shortening
- Excellent conformability and flexibility** to accommodate tortuous anatomy

### Covered stent

Indicated for the management of malignant and benign colorectal strictures



- Silicone coating** designed to resist tissue ingrowth and facilitate safe removal
- Retrieval string on distal end** to aid safe stent removal


## DESIGN CHOICE

The Gate Colorectal Stent is available in a broad range of sizes, with either a short or long delivery system. Please contact us with your preferred stent specification and we will provide you with additional ordering information.


STENT TYPE	DELIVERY SYSTEM		STENT	
	Diameter (Fr)	Length (mm)	Diameter (mm)	Length (cm)
UNCOVERED / STRAIGHT	10	700, 2200	18, 20, 22, 24, 26	4, 5, 6, 7, 8, 9, 10, 11, 12
COVERED / FLARED	10	700, 2200	18, 20	4, 5, 6, 7, 8, 9, 10, 11, 12

For additional information, please call us on [01844 340 620](tel:01844340620) or email [info@aprmedtech.com](mailto:info@aprmedtech.com)


Lumen conformity



Shape flexibility



Radial strength



Manufactured by:



[www.mediwood.co.kr](http://www.mediwood.co.kr)

Supplied exclusively in the UK by:



[aprmedtech.com](http://aprmedtech.com)  
@aprmedtech

APR Medtech Limited  
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)

