


# Apollo™

## Non-Compliant PTCA Balloon Catheter


A non-compliant high pressure balloon catheter designed for consistent performance and reliability

- Pole Vault™ shaft technology for optimal pushability and flexibility
- Laser cut, spiral hypotube for excellent trackability
- Well balanced tip for zero transition
- Minimal 3% axial balloon growth and low 0.55% compliance
- 22 atm Rated Burst Pressure
- Slider™ hydrophilic external surface from tip to adjacent guidewire port



### TECHNICAL SPECIFICATION

Balloon diameter (mm)	2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0, 4.5, 5.0
Balloon length (mm)	8, 10, 12, 15, 18
Nominal pressure	12 atm
Rated burst pressure	22 atm (Ø 2.0 - 4.0) 20 atm (Ø 4.5- 5.0)
Average burst pressure	30 atm
Compliance	Non-compliant, 0.55%
Balloon pleats	3
Tip length	2.0 / 2.5mm
Distal shaft OD	2.55F / 2.6F
Proximal shaft OD	2.0F
Working length (cm)	140cm
Guidewire compatibility	0.014"
Guiding catheter compatibility	5F compatible / 6F compatible for kissing
Coating	Slider™ hydrophilic external surface from tip to adjacent guidewire port



## ORDERING INFORMATION

For additional information, or to order one of our products, please call us on 01844 340 620 or email [info@aprmedtech.com](mailto:info@aprmedtech.com).

Diameter (mm)	Length				
	8	10	12	15	18
2.0	802-2008	802-2010	802-2012	802-2015	802-2018
2.25	802-2208	802-2210	802-2212	802-2215	802-2218
2.5	802-2508	802-2510	802-2512	802-2515	802-2518
2.75	802-2708	802-2710	802-2712	802-2715	802-2718
3.0	802-3008	802-3010	802-3012	802-3015	802-3018
3.25	802-3208	802-3210	802-3212	802-3215	802-3218
3.5	802-3508	802-3510	802-3512	802-3515	802-3518
3.75	802-3708	802-3710	802-3712	802-3715	802-3718
4.0	802-4008	802-4010	802-4012	802-4015	802-4018
4.5	802-4508	802-4510	802-4512	802-4515	802-4518
5.0	802-5008	802-5010	802-5012	802-5015	802-5018

Manufactured by:



[www.brosmed.com](http://www.brosmed.com)

Distributed in the United Kingdom by:



[aprmedtech.com](http://aprmedtech.com)  
@aprmedtech

APR Medtech Limited  
The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)

