


# ARTIMES™

## Semi-Compliant PTCA Balloon Catheter

A semi-compliant 'workhorse' balloon catheter designed to offer enhanced flexibility and crossability

- Ultra low crossing profile: 0.0186"
- Well balanced tip for zero transition
- Broad range of diameters and lengths
- Slider™ hydrophilic external surface from tip to adjacent guidewire port
- Eel™ hydrophobic coating inside guidewire lumen



### TECHNICAL SPECIFICATION

Balloon diameter (mm)	1.0, 1.25, 1.5, 1.75, 2.0, 2.25, 2.5, 2.75, 3.0, 3.25, 3.5, 3.75, 4.0
Balloon length (mm)	5, 8, 10, 12, 15, 20, 25, 30
Nominal pressure	6 atm
Rated burst pressure	14 atm
Average burst pressure	20 atm
Compliance	Semi-compliant
Balloon pleats	2/3
Tip length	1.5 / 2.0 / 2.5mm
Distal shaft OD	2.36F / 2.55F / 2.7F
Proximal shaft OD	1.9F
Working length	140cm
Guidewire compatibility	0.014"
Guiding catheter compatibility	5F compatible / 6F compatible for kissing
Coating	Slider™ hydrophilic external surface from tip to adjacent guidewire port Eel™ hydrophobic coating inside guidewire lumen



## ORDERING INFORMATION

For additional information, or to order one of our products, please call us on 01844 340 620 or email [info@aprmedtech.com](mailto:info@aprmedtech.com).

Diameter (mm)	Length							
	5	8	10	12	15	20	25	30
1.0	801-1005	801-1008	801-1010	801-1012	801-1015	801-1020	801-1025	801-1030
1.25	801-1205	801-1208	801-1210	801-1212	801-1215	801-1220	801-1225	801-1230
1.5	801-1505	801-1508	801-1510	801-1512	801-1515	801-1520	801-1525	801-1530
1.75	801-1705	801-1708	801-1710	801-1712	801-1715	801-1720	801-1725	801-1730
2.0	801-2005	801-2008	801-2010	801-2012	801-2015	801-2020	801-2025	801-2030
2.25	801-2205	801-2208	801-2210	801-2212	801-2215	801-2220	801-2225	801-2230
2.5	801-2505	801-2508	801-2510	801-2512	801-2515	801-2520	801-2525	801-2530
2.75	801-2705	801-2708	801-2710	801-2712	801-2715	801-2720	801-2725	801-2730
3.0	801-3005	801-3008	801-3010	801-3012	801-3015	801-3020	801-3025	801-3030
3.25	801-3205	801-3208	801-3210	801-3212	801-3215	801-3220	801-3225	801-3230
3.5	801-3505	801-3508	801-3510	801-3512	801-3515	801-3520	801-3525	801-3530
3.75	801-3705	801-3708	801-3710	801-3712	801-3715	801-3720	801-3725	801-3730
4.0	801-4005	801-4008	801-4010	801-4012	801-4015	801-4020	801-4025	801-4030

Manufactured by:



[www.brosmed.com](http://www.brosmed.com)

Distributed in the United Kingdom by:



[aprmedtech.com](http://aprmedtech.com)  
[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited  
 The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX  
 t. +44 1844 340 620 f. +44 1844 340 621 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)

