

# BARRIER MEDICAL

## Unique design. High protection.


The Barrier Medical Armboard has been carefully designed to cover the gaps that other boards miss and provide excellent operator protection from harmful scatter radiation.

### Patient arm support

Shielding moved inwards to provide unobstructed access for the procedure


### Mattress anchor

placed under the patient and mattress to secure the device



### Lead lining


reduces radiation scatter to the operating team by up to 96%\*



## ORDERING INFORMATION

For additional information, or to request an on-site demo of this device, please call us on 01844 340 620 or email [info@aprmedtech.com](mailto:info@aprmedtech.com)

Description	Product code	Units per box
Barrier Medical Armboard	BARR1	1



Manufactured by:



[www.barriermedical.co.uk](http://www.barriermedical.co.uk)

Supplied exclusively in the United Kingdom by:



[aprmedtech.com](http://aprmedtech.com)

[@aprmedtech](https://twitter.com/aprmedtech)

APR Medtech Limited

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX

t. +44 1844 340 620 e. [info@aprmedtech.com](mailto:info@aprmedtech.com)

Scan to see the data

