




Haemorrhage Occluder Pin with applicator



'Presacral bleeding during pelvic surgery can be a very severe problem. When it cannot be controlled by conventional methods, it can result in massive hemorrhage and even death.'¹

Haemorrhage Occluder Pin


A simple solution to a very severe problem

Serrated titanium

A biocompatible material commonly used for internal implants.² The serrated shaft of the pin offers a secure placement into the sacrum.

Choice of sizes

Available in 10mm & 14mm to tackle bleeds of varying size and location.



Malleable applicator for efficient placement

The Haemorrhage Occluder Pin comes preloaded in an easy to use applicator. The long shaft of the applicator is malleable and provides optimal visibility for placement.

Manufactured by:



SUR-HOP-COLO-A5-v2-1125.

Ordering information

For additional information, or to order this product, please call our Customer Support Team on **01844 340 620** or email info@aprmtech.com

Code	Description	Box qty
CR1007	Haemorrhage Occluder Pin Standard - 10mm	2
CR1014	Haemorrhage Occluder Pin Large - 14mm	2

Additional product information is available at:

aprmtech.com

1. Wang Q., Shi W., Zhou W., He Z., "New Concepts in Severe Presacral Haemorrhage During Proctectomy," Arch Surg., 1985;1120:1013-1020.
2. Dougherty, S.H., Simmon R.L.: "Infections in Bionic Man: The Pathobiology of Infections in Prosthetic Devices," Curr. Probl. Surg., 1982: 221-264

aprmtech



aprmtech.com

APR Medtech Ltd

The Sanderum Centre, 30a Upper High Street, Thame, Oxfordshire, OX9 3EX

t. +44 (0)1844 340 620

f. +44 (0)1844 340 621

e. info@aprmtech.com